


Ablation Pattern

Fractal Ablation

Refractive ablation






Full optical zone.
Percentage correction.

Small optical zone. 100% refractive.

Beam Pattern

Large angle Scanning

Parallel Scanning




Laser eyepiece



- Close-Up surgery.
- 3 cm working distance only.
- Stable and well controlled working distance.
- Integrated plume evacuation.
- Perfect alignment with visual axis.

Fully automated procedure



- Simplifies LASIK / PRK procedure down to the push of a button.
- Automated pupil detection.
- Video eyetracking system.

Video-eyetracker



Pure software eye tracking at a 1 mm tracking range and larger auto centration range

Further details on the M7

- Random projection for the best arrangement of laser shots on time and spatial domain.
- Fully automated arm for laser delivery that allows plenty of space to prepare the eye.
- Non-movable bed for perfect comfort of the patient.
- Small footprint of the Laser as a whole (including bed).